

# Penlon Sigma Delta Vaporizer

ANAESTHESIA SOLUTIONS

- ◆ Multiple agent and filler block options
- ◆ Five backbar mounting systems, including Dräger Auto-Exclusion®†
- ◆ Superb performance, particularly at low flow
- ◆ Service Free\*
- ◆ Low body weight



# Penlon Sigma Delta Vaporizer

ANAESTHESIA SOLUTIONS




The Sigma Delta vaporizer delivers accurate concentrations under varying conditions of flow rate and temperature, particularly at low flow


Vaporizer Type	Sevoflurane 8% Filler Type			Isoflurane 5% Filler Type		Halothane 5% Filler Type		Enflurane 5%, 7% Filler Type	
	Pour fill	Key fill	Quik-Fil	Pour fill	Key fill	Pour fill	Key fill	Pour fill	Key fill
	Selectatec Compatible with Interlock	✓	✓	✓	✓	✓	✓	✓	✓
Dräger Plug-In Interlock Compatible	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dräger Auto-Exclusion® Compatible*	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cagemount Taper 23 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓
North American Dräger Compatible	✓	✓	✓	✓	✓	✓	✓	✓	✓

- ◆ Key fill, Quik Fil®, or Pour fill
- ◆ Sevoflurane, Isoflurane, Halothane and Enflurane
- ◆ Selectatec®, Dräger Plug-In®, Dräger Auto-Exclusion®\*, North American Dräger, Cagemount
- ◆ Internationally recognised, award-winning design

## Typical performance graph



## Backbar mounting systems



Like us at [facebook.com/PenlonGlobal](https://facebook.com/PenlonGlobal)



Follow us at [twitter.com/PenlonGlobal](https://twitter.com/PenlonGlobal)

Doc. No. SD02175L / December 2017. Penlon, InterMed and Sigma Delta are trademarks of Penlon Limited. All other trademarks are the property of their respective owners. \*Presumes ten year product life requiring no preventive maintenance service. It is recommended that a service is carried out at ten years. Halothane vaporizers require a five year service. \*Auto-Exclusion not available in the US.

Part of the InterMed Group © Penlon Ltd. 2017 All rights reserved

Penlon Limited  
Abingdon Science Park  
Barton Lane, Abingdon  
OX14 3NB, UK

**General**  
t +44 (0) 1235 547000  
f +44 (0) 1235 547041  
w [www.penlon.com](http://www.penlon.com)

**International Sales**  
t +44 (0) 1235 547001  
f +44 (0) 1235 547021  
e [international.sales@penlon.com](mailto:international.sales@penlon.com)

**UK Sales**  
t +44 (0) 1235 547036  
f +44 (0) 1235 547023  
e [uk.sales@penlon.com](mailto:uk.sales@penlon.com)

**Technical Support**  
t +44 (0) 1235 547060  
f +44 (0) 1235 547061  
e [tech.support@penlon.com](mailto:tech.support@penlon.com)

