

# Penlon Sigma Delta Vaporizer

ANESTHESIA SOLUTIONS

- ◆ Multiple agent and filler block options
- ◆ Five backbar mounting systems, including Dräger Auto-Exclusion<sup>®†</sup>
- ◆ Superb performance, particularly at low flow
- ◆ Service Free\*
- ◆ Low body weight



# Penlon Sigma Delta Vaporizer


ANESTHESIA SOLUTIONS



The Sigma Delta vaporizer delivers accurate concentrations under varying conditions of flow rate and temperature, particularly at low flow


Vaporizer Type	Sevoflurane 8%			Isoflurane 5%		Halothane 5%		Enflurane 5%, 7%	
	Filler Type			Filler Type		Filler Type		Filler Type	
	Pour fill	Key fill	Quik-Fil	Pour fill	Key fill	Pour fill	Key fill	Pour fill	Key fill
Selectatec Compatible with Interlock	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dräger Plug-In Interlock Compatible	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dräger Auto-Exclusion <sup>+</sup> Compatible	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cagemount Taper 23 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓
North American Dräger Compatible	✓	✓	✓	✓	✓	✓	✓	✓	✓

## Typical performance graph



- ◆ Key fill or Pour fill
- ◆ Sevoflurane, Isoflurane, Halothane and Enflurane
- ◆ Selectatec<sup>®</sup>, Dräger Plug-In<sup>®</sup>, Dräger Auto-Exclusion<sup>®+</sup>, North American Dräger, Cagemount
- ◆ Internationally recognized, award-winning design

## Backbar mounting systems



Like us at [facebook.com/PenlonGlobal](https://facebook.com/PenlonGlobal)



Follow us at [twitter.com/PenlonGlobal](https://twitter.com/PenlonGlobal)

Doc. No. SD0217SL-US / December 2017. Penlon, InterMed and Sigma Delta are trademarks of Penlon Limited. All other trademarks are the property of their respective owners. \*Presumes ten year product life requiring no preventive maintenance service. It is recommended that a service is carried out at ten years. Halothane vaporizers require a five year service. <sup>+</sup>Auto-Exclusion not available in the US.

Part of the InterMed Group © Penlon, Inc. 2017 All rights reserved

Penlon, Inc.  
11515 K-Tel Drive  
Minnetonka  
MN 55343

**General**  
phone 952.933.3940  
fax 952.933.3375  
www.penlon.com

**US Sales**  
toll free 800.328.6216  
toll free fax 888.223.9540  
us.sales@penlon.com

**US Service**  
fax 952.697.2038  
us.service@penlon.com  
customer.service@penlon.com

**International Sales**  
phone 011 44 1235 547 001  
fax 011 44 1235 547 021  
international.sales@penlon.com

